

Highest precision without grid connection

Powerful and mobile

The rechargeable handpiece Q Mobile combines two worlds for the first time: unlimited mobility thanks to battery power and the quality of a Schick handpiece.

Your benefits

- ► Highest precision and concentricity properties
- ▶ High performance of 5.5 Ncm and 30,000 rpm
- Quick-release clamping
- Speed and capacity display
- ▶ Left- and right hand motion
- Optional: connecting cable for separation of battery and handpiece





Q Mobile motorhandpiece without cable, control unit with integrated battery, charging station, transport case

Technical Data	
Speed range	3,000-30,000 rpm
Torque	5.5 Ncm
Operating voltage	90-264 Volt
Operating frequency	50/60 Hz
Overload protection	Electrical

Handpiece	
Diameter	Max. 29 mm
Length without cable	120 mm
Weight without cable	194 g
Drive system	Brushless DC-Motor
Concentricity	< 0,02 mm
Chucks	Ø 2.35 mm (Standard) Ø 3.0 mm (Optional)
Tool change	Ouick release

Battery	
Width	71 mm
Height	125 mm
Depth	43 mm
Weight	250 g
Capacity	3,500 mAh



Mobile Sets



Q Mobile + connecting cable

Q Mobile complete Art.-No. 9705 with connecting cable Art.-No. 9604

Scope of delivery: Q Mobile motorhandpiece, control unit with integrated battery, charging station, connecting cable, transport case



Battery set for Q Motor with ISO-coupling

Art.-No. 9705/1

Scope of delivery: Q Motor with ISO-coupling, contol unit with integrated battery, charging station, connecting cable, transport case (contra-angle handpiece not included)



Battery set for your Schick motorhandpiece

Art.-No. 9988

Suitable for following motorhandpieces: Q Basic and Slim (Art.-No. 9400; 9400/3; 9400/10 and 9400/13), Q Profi (Art.-No. 9310 and 9300), QUBE LV (Art.-No. 9002), C2 Master (Art.-No. 6600), Q motor with ISO-coupling (Art.-No. 9200/04)

Scope of delivery: control unit with integrated battery, charging station, connecting cable, transport case

Optionally available:



► connecting cable Art.-Nr. 9604 (for Q Mobile) Art.-Nr. 9604/1 (for Q Basic, Q Profi, QUBE LV, C2 Master, Q Motor with ISO-coupling)



Schick GmbH

Lehenkreuzweg 12 D-88433 Schemmerhofen Telefon +49 7356 9500-0 Telefax +49 7356 9500-95 info@schick-dental.de www.schick-dental.de